


2007

Important

Before printing, turn off "fit to page" in your print options dialogue box.*




Length Measurement

- 1**
 - Place bottom edge of this sheet against wall on a hard surface.
 - Tape down to prevent any mis-sizing
- 2**
 - Place heel of AFO on sizing chart **against wall**
 - Be certain to include posterior stop
 - If foot is longer than brace, include foot
- 3**
 - Trace AFOs for accuracy
 - Match longest point to number on right
 - If you are between or slightly above line, round up

4 IMPORTANT
Necessary Adjustments for Circumference

For all AFOs, TCFlex, DAFOs, KAFOs, SMOs with:

- Solid ankle, instep strapping
Increase 1 full size
- Articulating hinges, posterior stop and instep strapping
Increase 1.5
- Turbo boot, posterior stop, hinges, footplates, instep, forefoot and toe strapping:
Increase 2 sizes



PLACE HEEL OF AFO AGAINST WALL



children and adolescent sizes



Keeping Pace, Inc.
contact us:
 footwear@keepingpace.com
 (888) 526-0020
 14 Fuller St.
 Gloucester, MA 01930
 USA
visit us on the web:
 www.keepingpace.com